


This note explains various steps involved in de-convoluting the right target(s) of small molecule by immobilizing the bait-molecule on solid support.

### Step-1) Derivatization of 'bait-molecule'

- 1) Analyze SAR (Structure-Activity Relationship) of the bait molecule.
- 2) Define a point of derivatization for minimal activity loss
- 3) Synthesis a primary amine (or other reactive functional group) derivative of the compound



### Step-2) Immobilization of 'bait-molecule' on Agarose Beads

'Bait-molecule' specific affinity matrix is prepared by immobilizing primary amine derivative of the compounds on epoxy activated agarose beads. Epoxy group on the beads are placed on a 12 atom long hydrophilic spacer for efficient coupling of small-molecule ligands. 'Control' affinity matrix is prepared by blocking the activated epoxy group through amine-based buffers.

### Step-3) Affinity Chromatography and Isolation of Target Proteins

- 1) Well washed 'Bait-molecule' specific and control affinity matrix are incubated with cell-lysate (protein source) for 3 hours in high salt buffer (1M NaCl in PBS) at 4 degree centigrade.
- 2) After incubation beads are briefly centrifuged and unbound protein fraction (supernatant) is separated.
- 3) Matrix is then washed 2X times with high salt buffer and 2X times with low salt buffer (150 mM NaCl in PBS).
- 4) Proteins bound to beads are then eluted 2X with elution buffer (1mM bait-molecule in PBS).
- 5) Proteins from the elution are then precipitated by chloroform-methanol based method.

### Step-4) Protein Identification

- 1) Precipitated proteins are separated over SDS-PAGE gel.
- 2) Separated proteins are then subjected to 'in-gel' trypsin digestion protocol and proteins are identified using mass-spectrometry based workflow.

### Step-5) Deconvolution of Targets and Conclusion

Target capture experiments are performed in triplicates. Specific target proteins are de-convoluted by comparing the protein profile obtained from bait-molecule specific affinity matrix and control experiments. Proteins that are significantly enriched in test experiments compared to control experiments are considered as specific binding partners of the test molecule.

## Typical Target Deconvolution Approach in 'Target Enrichment' Based Chemical-Proteomics

