


# Various Target Identification Technologies: Considerations and Eventual Deliverables

Disclaimers:

- 1) This table provides our general understanding of different technologies as learnt from publicly available literature.
- 2) Shantani do not practice or claim to practice all the methods of target identifications that are mentioned here.
- 3) UPT and SCLS are Shantani's Proprietary Technologies.



## Abbreviation

**UPT** = Unique Polymer Technology  
(Shantani's Proprietary Technology)

**TPP** = Thermal Proteome Profiling  
*Savitski et al., Science (2014) 346 (6205):1225784*

**DARTS** = Drug Affinity Responsive Target Stability  
*Lomenick et al. PNAS (2009) 106:21984-9*

**ABPP** = Activity Based Protein Profiling  
*Evans et al. Nat Biotechnol (2005) 23:1303-7*

**SCLS** = Subcellular Location Specific Target Capture & Identification Technology  
(Shantani's Proprietary Technology)

For Other Available Target Identification Methodologies:

*Saxena et al. Expert Opin. Drug Discov. (2009) 4:701-714*

*Saxena C. Expert Opin. Drug Discov. (2016) 11:1017-1025*



# Target Identification: Comparison of Chemical-Proteomics Based Mass-Spectrometry Based Methodologies

	Resources	Molecule Derivatization Not Needed / Structure-Activity Relationship not Available				Molecule Needs to be Derivatized / Structure-Activity Relationship Available			
Criteria		Computational Tools	UPT	TPP	DARTS	Bead / Biotinylation	Photo-Cross Linking	ABPP	SCLS
Resources	Chemistry	-	-	-	-	++	+++	+++	++
	Computational and/OR Statistical	+++++	++	++++	++++	++	+++	+++	++
	Biochemistry	-	+	+	+	+	+	+	++
	Cell-Biology	-	+	+	+	+	+	+	++
	Mass-Spectrometry	-	++	+++	+++	+++	+++	+++	++
Time		+	++	++	++	++++	++++	++++	+++++
Cost		+	+++	++++	++++	+++	++++	++++	+++++
False Positive Identification		~ 80 %	~40%	~60%	~60%	~40%	~40%	~30%	~20%
Key Advantage(s)		Quick & Cost Effective	<ul style="list-style-type: none"> <li>Underivatized bait-molecule</li> <li>Target Identification through affinity based target enrichment</li> </ul>	Underivatized bait-molecule	Underivatized bait-molecule	Targets of varied affinity (high and low) are captured	Covalent capture of target	Activity based target capture	Capture of target from live cell in sub-cellular location specific manner
Disadvantages		Very High False Positive Rate	Certain highly water Soluble Compounds can not used	<ul style="list-style-type: none"> <li>Limited to Soluble Proteins only</li> <li>Complex data-analysis</li> </ul>	<ul style="list-style-type: none"> <li>Limited to Highly Abundant Proteins only</li> </ul>	<ul style="list-style-type: none"> <li>Bait-molecule derivatization is needed for immobilization</li> <li>Limited success in capturing membrane targets</li> </ul>	<ul style="list-style-type: none"> <li>Bait-molecule derivatization is needed for immobilization</li> <li>Primary bait-molecule may loose activity because relatively larger size of probes</li> </ul>	<ul style="list-style-type: none"> <li>Bait-molecule derivatization is needed for immobilization</li> <li>Very Specific probes needs to be synthesized –target class specific capture</li> </ul>	<ul style="list-style-type: none"> <li>Bait-molecule derivatization is needed for immobilization</li> <li>Particular Function / phenotype linked target are captured</li> </ul>

